

# Synchronizer

RF1-L

product manual

version 1.0

Dr.X Works BV  
Engelsbergenstraat 18  
5616JC Eindhoven  
The Netherlands

Phone: +31852460114  
[www.drx-works.nl](http://www.drx-works.nl)  
[info@drx-works.nl](mailto:info@drx-works.nl)



# Contents

<b>1 Introduction</b>	<b>4</b>
1.1 block diagram . . . . .	4
<b>2 Installation</b>	<b>7</b>
2.1 Unpacking and inspection . . . . .	7
2.2 Connections and controls . . . . .	7
2.3 power requirements . . . . .	7
<b>3 Operation</b>	<b>9</b>
3.1 Open/Closed loop . . . . .	9
3.2 Calibration procedure . . . . .	9
3.3 Daily operation . . . . .	10
3.4 Operation without laser source . . . . .	10
<b>4 Technical specification</b>	<b>11</b>
4.1 Electrical . . . . .	11
4.2 Photodiode detector . . . . .	11
4.3 Mechanical . . . . .	11
<b>5 Environmental</b>	<b>12</b>
<b>6 Maintenance</b>	<b>13</b>
<b>7 Safety</b>	<b>14</b>
<b>8 Warranty</b>	<b>15</b>
<b>9 Mechanical outline</b>	<b>16</b>
<b>Bibliography</b>	<b>17</b>
<b>A External phase control scheme</b>	<b>18</b>
<b>B Optional product: Synchronizer remote</b>	<b>19</b>
B.1 The package . . . . .	19
B.2 Specifications . . . . .	19
B.3 Programming . . . . .	19

# 1. Introduction

This document details the design of the Phase Locked Loop (PLL) Synchronizer. The PLL Synchronizer is designed to synchronize a 3-GHz RF oscillator to a mode-locked 75-MHz Ti:Sapphire laser system with less than 1-ps phase jitter. This PLL Synchronizer model is a commercial version of the first version that has been installed at Eindhoven University of Technology for the purpose of an RF photogun setup. The principle of the synchronization has been published in ref.[1].

In many experimental systems mode-locked lasers have been synchronized to other lasers, to electronic oscillators and, recently, to a free electron laser with subpicosecond jitter. In these applications the frequency of the laser is adjusted by changing the cavity length, to follow a master oscillator. In case the master oscillator is a stable crystal oscillator the advantages are obvious: the absolute (long term) jitter approaches the jitter of the crystal oscillator which is usually much better than the jitter and drift of a laser oscillator. However, in many applications it is of greater interest how well the two sources are synchronized, i.e. the relative jitter counts and not the absolute jitter. In pump-probe measurements, for instance, it is the delay between the pump and the probe which is of prime importance, irrespective of the absolute timing with respect to some external reference. In these cases changing the frequency of the oscillator with the highest control bandwidth is to be preferred, because in a stable well-designed feedback system this bandwidth will determine the residual jitter. The control bandwidth of electronic oscillators often exceeds the control bandwidth of laser oscillators by more than three orders of magnitude. Typically, the control bandwidth of the 75-MHz laser oscillator is of the order of 1 kHz compared to a few MHz for the 3-GHz electronic oscillator. By making the laser the master oscillator in the design of the PLL Synchronizer the relative phase-jitter is improved by almost 2 orders of magnitude, resulting in an RMS jitter of less than 50 fs.

## 1.1 block diagram

A schematic layout of the synchronization system is given in Fig. 1.1. The Ti:Sapphire laser system should operate at a repetition rate of 75 MHz. The repetition rate is inversely proportional to the cavity length and can therefore be changed by moving one of the end mirrors by for instance a piezo-electric transducer (PZT). The electronic oscillator of the PLL Synchronizer operates at a frequency of 3 GHz. Its frequency is voltage controllable with a bandwidth of 3 MHz. The 75-MHz optical pulse train from the Ti:Sapphire laser is monitored with a biased fast Si-(PIN)-photodiode. A narrow-band filter is used to select the 5th harmonic of this signal at 375 MHz. As phase-jitter increases at higher harmonics the use of the 5th harmonic increases the sensitivity of the loop. To compare the phase of

the laser with the phase of the electronic oscillator the 3-GHz RF signal is divided by 8 and subsequently filtered to remove the higher harmonics. The narrow-band filters are passive band pass devices having a bandwidth of 16 MHz. The filters behind the amplifiers are used to remove spurious amplifier output. To avoid phase errors by temperature changes the two signal paths to the inputs of the mixer are constructed fully identically. A Double-Balanced Mixer (DBM) is used to measure the phase difference of the two amplified signals. The DBM output signal is used by the loop electronics to change the frequency of the electronic oscillator.



## 2. Installation

### 2.1 Unpacking and inspection

1. After installation into 19"-rack or other, connect mains cable to the back panel.
2. Connect the photodiode to the input on the front panel. It is recommended to use semi-rigid coaxial cable (e.g. RG405).
3. Connect TTL signal to "ENABLE" input on the front panel.
4. It is optional to control remotely the phase setting via back panel connector using a 5k potentiometer or an external source. It is recommended to use 4-core shielded cable. (appendix B)
5. After switching on the PLL Synchronizer the light input on the photodiode should be adjusted to set the value on the ratio monitor to  $1.0 \pm 10\%$ .
6. Perform calibration procedure that is described in section 3.2.

### 2.2 Connections and controls

In the table below the input and output connectors on the front and back panel are summarized.

#### Front panel:

Photo Detector	Input	SMA	50 $\Omega$	+15V bias
RF Enable	Input	BNC	50 $\Omega$	TTL
3 GHz	Output	SMA	50 $\Omega$	10 dBm
Phase Error	Output	BNC	150 $\Omega$	+/- 3V, typical 0.1 mV RMS

#### Back panel:

Remote Phase Control	Input	DIN	5 pin 240° circular DIN plug (appendix A)
Mains power	Input	Euroconnector	110 or 230 V (auto switching)
Out of Range	Output	2 pins	relay contact, max 24 V and 50 mA

### 2.3 power requirements

The current source operates on single phase AC mains power from 100 V to 240V with a frequency between 47 and 63 Hz without any need for range selection. The primary fuses

are 4AT.



## 3. Operation

### 3.1 Open/Closed loop

The PLL Synchronizer can run in two modes: 1- fixed frequency (open), and 2- phase locked (closed). In the first case, the voltage controlled oscillator is driven by a constant voltage such that the output frequency is fixed to 2998.5 MHz. In the latter case the loop is closed and the phase (and frequency) is continuously adjusted to lock the phase. The mode can be selected by the switch on the front panel.

### 3.2 Calibration procedure

The calibration of the phase error output is described in this section. It is recommended to perform this procedure after installation and on a regular basis to check the performance of the PLL Synchronizer.

1. Connect the phase error output to an oscilloscope. (1M $\Omega$  input impedance)
2. Set the PLL Synchronizer in open loop.
3. Tune the external oscillator connected to the photodiode input to the central frequency of the PLL Synchronizer (375MHz). On the oscilloscope screen you will see the beat frequency between the two oscillators. The beat frequency should be in the order of one kHz.
4. Set the PLL Synchronizer in closed loop. Now you will see a noisy trace on the oscilloscope.
5. The RMS value of this signal, which is measured preferably with a digital multimeter with true-RMS capability, is the resulting phase error of the PLL Synchronizer.

The corresponding time error at 3 GHz is easily calculated by the following formula:

$$\tau_{rms} = \frac{U_{\text{phaseerror}}}{2\pi\alpha} \frac{8}{f_{\text{RF}}} \quad (3.1)$$

6. The calibration factor ( $\alpha$ ) is expressed in V/rad. The phase difference between the positive and negative top is exactly  $\pi$  rad. Note that this value is the phase error calibration at 375 MHz.

### 3.3 Daily operation

If RF is disabled the RF power is dumped in the internal load of the RF switch. In case the RF is enabled, make sure that the 3-GHz RF output is always connected to a  $50\Omega$  load. For the stability of operation the signal strength should be checked regularly and if necessary the light intensity on the photodiode should be adjusted. Keeping the room temperature stable improves the performance of the PLL Synchronizer.

### 3.4 Operation without laser source

It is possible to operate the PLL Synchronizer without the mode-locked laser source, while still controlling the frequency of the PLL Synchronizer. In this case the photodiode cable has to be disconnected; instead a sine-wave generator of 375 MHz (+/- 0dBm) should be connected using a decoupling capacitor. In this case, the sine-wave generator acts as the master oscillator to which the electronic oscillator is phase locked.

## 4. Technical specification

### 4.1 Electrical

The PLL Synchronizer is mounted in an EMC enclosure.

Central Frequency	2998.5 MHz
Frequency range	+/- 1.5 MHz
Output power	10 dBm typical
RF enable	BNC connector, TTL
Phase control, Local	Dial, minimal $0 - 2\pi$ rad
Phase control, Remote	5 pin 240° circular DIN plug (appendix A), minimal $0 - 2\pi$ rad
Ratio Monitor	Digital front panel meter
Switch loop open/close	Constant frequency (2998.5 MHz) in case opened
In-range relay	Back panel connectors, max. +24 V and 50 mA
Photo detector bias	SMA connector, +15 V
Power Supply	Mains 80-130V for 60 Hz or 180-250V for 50 Hz

### 4.2 Photodiode detector

The PLL Synchronizer is supplied with an ultrafast silicon photodiode detector: type no. AEPX65 (manufacturer Centronic). The diode is biased with +15 V supplied by the PLL Synchronizer through the SMA connector.

### 4.3 Mechanical

Dimensions (LHW)	19" x 3U x 500mm
Weight	ca. 9kg

If installed into a 19" rack/chassis the PLL Synchronizer must be supported by its rear side. The front panel mounts are not sufficient alone to bear the weight of the PLL Synchronizer.

## 5. Environmental

Operating temperature range 0 to 40°C ambient, however temperature fluctuations should be kept as small as possible to prevent drifting. The Synchronizer should be operated in a clean environment. On no account must the Synchronizer be exposed to water ingress.

## 6. Maintenance

No maintenance is required in normal operation and there are no user-serviceable parts within the Synchronizer. Each product is factory tested and supplied with a set of test results. If degradation in performance to below the specified levels occurs, or a failure is suspected, then the complete unit should be returned to the manufacturer together with details of the fault.

## 7. Safety


Hazardous voltages are present within this line-operated unit. Do not remove any panels. The Synchronizer must be grounded. Line supply is filtered within the amplifier. A replaceable fuse is fitted to the line input connector which should be changed only when the AC power is disconnected from the connector.

## 8. Warranty


Dr.X Works B.V. warrants for 2 years from date of shipment that the goods supplied will be in full compliance with the agreed specifications and will be free from defects in material and workmanship. Any and all other warranties (except of title) express or implied, relating to fitness for particular purpose, merchantable quality or otherwise are expressly disclaimed. Seller will not be responsible for special or consequential loss or damages. Liability shall be limited to the repair or replacement of defective products subject to the return of the product intact, and un-tampered with by the buyer.

# 9. Mechanical outline

## Front panel



## Rear panel







# Bibliography

- [1] F. B. Kiewiet et al. “Femtosecond synchronization of a 3 GHz RF oscillator to a mode-locked Ti:sapphire laser”. In: *Nucl. Instrum. Meth. A* 484 (2002), p. 619. DOI: 10.1016/S0168-9002(01)01994-5.

## A. External phase control scheme



The shield of the cable is connected to pin 3 of the 5 pin 240° DIN-plug. Recommended potentiometer type: Spectrol/Vishay 534-1-1-5k 5%

## B. Optional product: Synchronizer remote

### B.1 The package

1. Remote extension lead (5 way, shielded) (to be delivered in user specified length)
2. USB cable
3. Synchronizer Remote control

### B.2 Specifications

Size (lxbxh)	169 x 165 x 53 mm
Weight	??

Typical performance, measured fully warmed-up (30 min.)

	usb	potentiometer
Full scale error(1)	-0.7%	+/- 5%
Resolution	12 bit	4000
Nonlinearity(1)	< 0.2%	< 0.25%
Temperature drift(1)	8.5 ppm/°C	20 ppm/°C

Notes: (1) these values represent the contribution of the remote control to the whole system.

### B.3 Programming

The address of the DAC is : 0x18 The busnumber of the DAC is 0x02 De virtual serial port has to be configured as follows:

Baudrate	57600
Parity	none
Databits	8
Stopbits	1
Hardware flow control	RTS/CTS

When the remote control is going to be used in USB mode, a virtual serial port has to be installed: <http://www.ftdichip.com/>. The windows device manager will show what portnumber is assigned to the remote control. To test the configuration free software is available at Sourceforge.com <http://sourceforge.net/projects/realterm/files/>.

To set a value to the DAC (eg. FFF) use the following commands:

Ascii	
"G"	Set busnummer
"2"	to "1"
"S"	Start I2C = "S" */
"1"	1st nibble adres
"8"	2nd nibble adres
"0"	1th nibble data (has to be 0)
"F"	2th nibble data
"F"	3th nibble data
"F"	4th nibble data
"P"	stop I2C